



Coagulation Testing

Reference Range Update

Client Services

Oahu808.677.7998

Neighbor
Islands.....866.281.6816

Effective November 20, 2018, Clinical Labs of Hawaii will be adjusting the Fibrinogen and Heparin Therapeutic Range. The PT and PTT reference ranges are unchanged.

The table below shows the previous and newly updated ranges. Prior to a change in reagent lot numbers, Clinical Labs of Hawaii performs extensive in-house studies to validate reagent lot performance, evaluate reagent lot consistency, and verify current ranges. Based on the results of these studies, it may be necessary at times to adjust our reference and therapeutic ranges.

Test	Current Range	New Range
Fibrinogen	205-445 mg/dL	214-493 mg/dL
Heparin Therapeutic Range	91-125 Seconds	81-137 Seconds

If you have any questions, please contact our Client Services Department.

Thank you for choosing Clinical Labs of Hawaii.